

GMIA STANDARD METHODS FOR THE TESTING OF EDIBLE GELATIN



2.7 Sulfur Dioxide

PRINCIPLE

Sulfur dioxide content can be determined using the Monier-Williams test method.

REFERENCES

AOAC Method 962.16 as stated in Food Chemicals Codex, General Tests and Assays, Appendix X.

REAGENTS AND SOLUTIONS

Stated in Sulfur Dioxide Determination in Food Chemicals Codex, General Tests and Assays, Appendix X

APPARATUS

Monier-Williams apparatus as per AOAC Method 962.16

PROCEDURE

Performed by contact laboratories as per AOAC Method 962.16.